

Product datasheet for **GA102577**

Human KIR2DS1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	KIR2DS1 CRISPRa kit - CRISPR gene activation of human killer cell immunoglobulin like receptor, two Ig domains and short cytoplasmic tail 1
Symbol:	KIR2DS1
UniProt ID:	Q14954
Synonyms:	CD158a; CD158H; p50.1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)