

Product datasheet for **GA102579**

Human KIR2DS4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	KIR2DS4 CRISPRa kit - CRISPR gene activation of human killer cell immunoglobulin like receptor, two Ig domains and short cytoplasmic tail 4
Symbol:	KIR2DS4
UniProt ID:	P43632
Synonyms:	CD158I; KIR1D; KIR2DS1; KIR412; KKA3; NKAT-8; NKAT8
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)