


## Product datasheet for **GA102583**

### Human KIR3DS1 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	KIR3DS1 CRISPRa kit - CRISPR gene activation of human killer cell immunoglobulin like receptor, three Ig domains and short cytoplasmic tail 1
Symbol:	KIR3DS1
Components:	GA102583G1, KIR3DS1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA102583G2, KIR3DS1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA102583G3, KIR3DS1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q14943</u>
Synonyms:	CD158E2; KIR-123FM; KIR-G1; NKAT-10; NKAT10
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online »](#)