

Product datasheet for **GA102603**

Human KPNA3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	KPNA3 CRISPRa kit - CRISPR gene activation of human karyopherin subunit alpha 3
Symbol:	KPNA3
UniProt ID:	<u>O00505</u>
Synonyms:	hSRP1; IPOA4; SRP1; SRP1gamma; SRP4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)