

Product datasheet for **GA102645**

Human Lymphocyte Activation Gene 3 (LAG3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	LAG3 CRISPRa kit - CRISPR gene activation of human lymphocyte activating 3
Symbol:	LAG3
UniProt ID:	P18627
Synonyms:	CD223
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)