

Product datasheet for **GA102650**

Human LAMA3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	LAMA3 CRISPRa kit - CRISPR gene activation of human laminin subunit alpha 3
Symbol:	LAMA3
UniProt ID:	Q16787
Synonyms:	BM600; E170; LAMNA; LOCS
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)