

## Product datasheet for **GA102672**

### Human LDHA activation kit by CRISPRa

#### Product data:

|               |  |
|---------------|--|
| Product Type: | CRISPRa Kits   |
| Description:  | LDHA CRISPRa kit - CRISPR gene activation of human lactate dehydrogenase A |
| Symbol:       | LDHA   |
| UniProt ID:   | <u><a href="#">P00338</a></u>  |
| Synonyms:     | GSD11; HEL-S-133P; LDHM; PIG19   |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)    |

#### Product images:



[View online >](#)