

Product datasheet for **GA102673**

Human Lactate Dehydrogenase B (LDHB) activation kit by CRISPRa

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | LDHB CRISPRa kit - CRISPR gene activation of human lactate dehydrogenase B |
| Symbol: | LDHB |
| UniProt ID: | P07195 |
| Synonyms: | HEL-S-281; LDH-B; LDH-H; LDHBD; TRG-5 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)