

Product datasheet for GA102675

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Human LDL Receptor (LDLR) activation kit by CRISPRa

Product data:

Product Type: CRISPRa Kits

Description: LDLR CRISPRa kit - CRISPR gene activation of human low density lipoprotein receptor

Symbol: LDLR UniProt ID: P01130

Synonyms: FH; FHC; LDLCQ2

Format: 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



