

## Product datasheet for **GA102677**

### Human Leptin (LEP) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	LEP CRISPRa kit - CRISPR gene activation of human leptin
Symbol:	LEP
UniProt ID:	<a href="#">P41159</a>
Synonyms:	LEPD; OB; OBS
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)