

Product datasheet for **GA102693**

Human LIFR activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	LIFR CRISPRa kit - CRISPR gene activation of human LIF receptor subunit alpha
Symbol:	LIFR
UniProt ID:	P42702
Synonyms:	CD118; LIF-R; SJS2; STWS; SWS
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)