

Product datasheet for **GA102695**

Human DNA Ligase III (LIG3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	LIG3 CRISPRa kit - CRISPR gene activation of human DNA ligase 3
Symbol:	LIG3
UniProt ID:	<u>P49916</u>
Synonyms:	LIG2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)