

Product datasheet for **GA102696**

Human DNA Ligase IV (LIG4) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	LIG4 CRISPRa kit - CRISPR gene activation of human DNA ligase 4
Symbol:	LIG4
UniProt ID:	<u>P49917</u>
Synonyms:	LIG4S
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)