

Product datasheet for **GA102734**

Human LRP3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	LRP3 CRISPRa kit - CRISPR gene activation of human LDL receptor related protein 3
Symbol:	LRP3
UniProt ID:	O75074
Synonyms:	low density lipoprotein receptor-related protein 3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)