

Product datasheet for **GA102738**

Human LRPAP1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	LRPAP1 CRISPRa kit - CRISPR gene activation of human LDL receptor related protein associated protein 1
Symbol:	LRPAP1
UniProt ID:	P30533
Synonyms:	A2MRAP; A2RAP; alpha-2-MRAP; HBP44; MRAP; MYP23; RAP
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)