

Product datasheet for **GA102739**

Human LSAMP activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	LSAMP CRISPRa kit - CRISPR gene activation of human limbic system associated membrane protein
Symbol:	LSAMP
Components:	GA102739G1, LSAMP gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA102739G2, LSAMP gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA102739G3, LSAMP gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q13449</u>
Synonyms:	IGLON3; LAMP
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)