

Product datasheet for **GA102742**

Human Leukotriene A4 hydrolase (LTA4H) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	LTA4H CRISPRa kit - CRISPR gene activation of human leukotriene A4 hydrolase
Symbol:	LTA4H
UniProt ID:	P09960
Synonyms:	leukotriene A4 hydrolase
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)