

Product datasheet for **GA102750**

Human LTC4S activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	LTC4S CRISPRa kit - CRISPR gene activation of human leukotriene C4 synthase
Symbol:	LTC4S
UniProt ID:	Q16873
Synonyms:	leukotriene C4 synthase; MGC33147
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)