

Product datasheet for **GA102836**

Human MC3R activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	MC3R CRISPRa kit - CRISPR gene activation of human melanocortin 3 receptor
Symbol:	MC3R
UniProt ID:	<u>P41968</u>
Synonyms:	BMIQ9; MC3; MC3-R; OB20; OQTL
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)