

Product datasheet for **GA102837**

Human MC4 Receptor (MC4R) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	MC4R CRISPRa kit - CRISPR gene activation of human melanocortin 4 receptor
Symbol:	MC4R
UniProt ID:	P32245
Synonyms:	MGC126851; MGC138197
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)