

Product datasheet for **GA102911**

Human Aquaporin 0 (MIP) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	MIP CRISPRa kit - CRISPR gene activation of human major intrinsic protein of lens fiber
Symbol:	MIP
UniProt ID:	P30301
Synonyms:	AQP0; CTRCT15; LIM1; MIP26; MP26
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)