

Product datasheet for **GA102922**

Human MLK3 (MAP3K11) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	MAP3K11 CRISPRa kit - CRISPR gene activation of human mitogen-activated protein kinase kinase kinase 11
Symbol:	MAP3K11
Components:	GA102922G1, MLK3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA102922G2, MLK3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA102922G3, MLK3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q16584
Synonyms:	MEKK11; MLK-3; MLK3; PTK1; SPRK
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)