

Product datasheet for **GA102924**

Human MLLT1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	MLLT1 CRISPRa kit - CRISPR gene activation of human MLLT1 super elongation complex subunit
Symbol:	MLLT1
Components:	GA102924G1, MLLT1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA102924G2, MLLT1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA102924G3, MLLT1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q03111</u>
Synonyms:	ENL; LTG19; YEATS1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)