

Product datasheet for **GA102935**

Human MMP3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	MMP3 CRISPRa kit - CRISPR gene activation of human matrix metallopeptidase 3
Symbol:	MMP3
UniProt ID:	P08254
Synonyms:	CHDS6; MMP-3; SL-1; STMY; STMY1; STR1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)