

Product datasheet for **GA102949**

Human MN1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	MN1 CRISPRa kit - CRISPR gene activation of human MN1 proto-oncogene, transcriptional regulator
Symbol:	MN1
Components:	GA102949G1, MN1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA102949G2, MN1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA102949G3, MN1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q10571
Synonyms:	dj353E16.2; MGCR; MGCR1; MGCR1-PEN
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)