

Product datasheet for **GA102976**

Human CITED1 activation kit by CRISPRa

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | CITED1 CRISPRa kit - CRISPR gene activation of human Cbp/p300 interacting transactivator with Glu/Asp rich carboxy-terminal domain 1 |
| Symbol: | CITED1 |
| UniProt ID: | Q99966 |
| Synonyms: | MSG1 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)