

Product datasheet for **GA102979**

Human MSH4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	MSH4 CRISPRa kit - CRISPR gene activation of human mutS homolog 4
Symbol:	MSH4
Components:	GA102979G1, MSH4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA102979G2, MSH4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA102979G3, MSH4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q15457</u>
Synonyms:	mutS homolog 4; mutS homolog 4 (E. coli); OTTHUMP00000011349
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)