

Product datasheet for **GA103003**

Human MT3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	MT3 CRISPRa kit - CRISPR gene activation of human metallothionein 3
Symbol:	MT3
UniProt ID:	P25713
Synonyms:	GIF; GIFB; GRIF; ZnMT3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)