

Product datasheet for **GA103034**

Human MUC3 (MUC3A) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	MUC3A CRISPRa kit - CRISPR gene activation of human mucin 3A, cell surface associated
Symbol:	MUC3A
Components:	GA103034G1, MUC3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA103034G2, MUC3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA103034G3, MUC3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q02505</u>
Synonyms:	MUC-3A; MUC3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)