

## Product datasheet for **GA103074**

### Human Myosin light chain 3 (MYL3) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	MYL3 CRISPRa kit - CRISPR gene activation of human myosin light chain 3
Symbol:	MYL3
UniProt ID:	<u><a href="#">P08590</a></u>
Synonyms:	CMH8; MLC1SB; MLC1V; VLC1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)