

Product datasheet for **GA103108**

Human NAP1L4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	NAP1L4 CRISPRa kit - CRISPR gene activation of human nucleosome assembly protein 1 like 4
Symbol:	NAP1L4
Components:	GA103108G1, NAP1L4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA103108G2, NAP1L4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA103108G3, NAP1L4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q99733
Synonyms:	hNAP2; NAP1L4b; NAP2; NAP2L
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)