

Product datasheet for **GA103190**

Human NFIL3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	NFIL3 CRISPRa kit - CRISPR gene activation of human nuclear factor, interleukin 3 regulated
Symbol:	NFIL3
UniProt ID:	Q16649
Synonyms:	E4BP4; IL3BP1; NF-IL3A; NFIL3A
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)