

Product datasheet for **GA103217**

Human Nkx3.1 (NKX3-1) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	NKX3-1 CRISPRa kit - CRISPR gene activation of human NK3 homeobox 1
Symbol:	NKX3-1
UniProt ID:	Q99801
Synonyms:	BAPX2; NKX3; NKX3.1; NKX3A
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)