


Product datasheet for **GA103236**

Human CCR4 NOT transcription complex subunit 3 (CNOT3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	CNOT3 CRISPRa kit - CRISPR gene activation of human CCR4-NOT transcription complex subunit 3
Symbol:	CNOT3
UniProt ID:	<u>Q75175</u>
Synonyms:	LENG2; NOT3; NOT3H
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)