

## Product datasheet for **GA103246**

### Human Nucleoside phosphorylase (PNP) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	PNP CRISPRa kit - CRISPR gene activation of human purine nucleoside phosphorylase
Symbol:	PNP
UniProt ID:	<a href="#">P00491</a>
Synonyms:	NP; PRO1837; PUNP
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)