

Product datasheet for **GA103260**

Human Natriuretic Peptide Receptor C (NPR3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	NPR3 CRISPRa kit - CRISPR gene activation of human natriuretic peptide receptor 3
Symbol:	NPR3
UniProt ID:	<u>P17342</u>
Synonyms:	ANP-C; ANPR-C; ANPRC; C5orf23; GUCY2B; NPR-C; NPRC
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)