

Product datasheet for **GA103300**

Human OAS3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	OAS3 CRISPRa kit - CRISPR gene activation of human 2'-5'-oligoadenylate synthetase 3
Symbol:	OAS3
UniProt ID:	Q9Y6K5
Synonyms:	p100; p100OAS
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)