

Product datasheet for **GA103321**

Human OPHN1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	OPHN1 CRISPRa kit - CRISPR gene activation of human oligophrenin 1
Symbol:	OPHN1
Components:	GA103321G1, OPHN1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA103321G2, OPHN1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA103321G3, OPHN1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	O60890
Synonyms:	ARHGAP41; MRX60; OPN1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)