

## Product datasheet for **GA103337**

### Human SLC22A18AS activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	SLC22A18AS CRISPRa kit - CRISPR gene activation of human solute carrier family 22 member 18 antisense
Symbol:	SLC22A18AS
Components:	GA103337G1, SLC22A18AS gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA103337G2, SLC22A18AS gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA103337G3, SLC22A18AS gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u><a href="#">Q8N1D0</a></u>
Synonyms:	BWR1B; BWSCR1B; ORCTL2S; p27-BWR1B; SLC22A1LS
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online »](#)