

Product datasheet for **GA103352**

Human P2X1 (P2RX1) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	P2RX1 CRISPRa kit - CRISPR gene activation of human purinergic receptor P2X 1
Symbol:	P2RX1
UniProt ID:	P51575
Synonyms:	P2X1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)