

## Product datasheet for **GA103353**

### Human P2X3 (P2RX3) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	P2RX3 CRISPRa kit - CRISPR gene activation of human purinergic receptor P2X 3
Symbol:	P2RX3
UniProt ID:	<a href="#">P56373</a>
Synonyms:	P2X3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)