

Product datasheet for **GA103354**

Human P2RX4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	P2RX4 CRISPRa kit - CRISPR gene activation of human purinergic receptor P2X 4
Symbol:	P2RX4
UniProt ID:	Q99571
Synonyms:	P2X4; P2X4R
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)