

Product datasheet for **GA103355**

Human P2X5 (P2RX5) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	P2RX5 CRISPRa kit - CRISPR gene activation of human purinergic receptor P2X 5
Symbol:	P2RX5
UniProt ID:	Q93086
Synonyms:	LRH-1; P2X5; P2X5R
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)