

Product datasheet for **GA103357**

Human P2Y1 (P2RY1) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	P2RY1 CRISPRa kit - CRISPR gene activation of human purinergic receptor P2Y1
Symbol:	P2RY1
UniProt ID:	P47900
Synonyms:	P2Y1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)