

Product datasheet for **GA103358**

Human P2Y2 (P2RY2) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	P2RY2 CRISPRa kit - CRISPR gene activation of human purinergic receptor P2Y2
Symbol:	P2RY2
UniProt ID:	P41231
Synonyms:	HP2U; P2RU1; P2U; P2U1; P2UR; P2Y2; P2Y2R
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)