

## Product datasheet for **GA103359**

### Human P2Y4 (P2RY4) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	P2RY4 CRISPRa kit - CRISPR gene activation of human pyrimidinergetic receptor P2Y4
Symbol:	P2RY4
Components:	GA103359G1, P2Y4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA103359G2, P2Y4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA103359G3, P2Y4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>P51582</u>
Synonyms:	NRU; P2P; P2Y4; UNR
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)