

Product datasheet for **GA103360**

Human P2Y6 (P2RY6) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	P2RY6 CRISPRa kit - CRISPR gene activation of human pyrimidinergetic receptor P2Y6
Symbol:	P2RY6
UniProt ID:	Q15077
Synonyms:	P2Y6
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)