

Product datasheet for **GA103361**

Human P2Y11 (P2RY11) activation kit by CRISPRa

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | P2RY11 CRISPRa kit - CRISPR gene activation of human purinergic receptor P2Y11 |
| Symbol: | P2RY11 |
| UniProt ID: | Q96G91 |
| Synonyms: | P2Y11 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)