

## Product datasheet for **GA103363**

### Human P4HB activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	P4HB CRISPRa kit - CRISPR gene activation of human prolyl 4-hydroxylase subunit beta
Symbol:	P4HB
UniProt ID:	<a href="#">P07237</a>
Synonyms:	DSI; ERBA2L; GIT; P4Hbeta; PDI; PDIA1; PHDB; PO4DB; PO4HB; PROHB
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)