

## Product datasheet for **GA103368**

### Human PACE4 (PCSK6) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	PCSK6 CRISPRa kit - CRISPR gene activation of human proprotein convertase subtilisin/kexin type 6
Symbol:	PCSK6
UniProt ID:	<a href="#">P29122</a>
Synonyms:	PACE4; SPC4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)